

One-Fil®

CALCIUM SILICATE-BASED BIOCERAMIC ROOT CANAL SEALER

FDA


Advantages

- High radiopacity
- High wettability
- High sealing ability
- High biocompatibility
- High dimensional stability

Setting time

- 30min in 100% humidity condition. (ISO6876:2012)
- 3.5 hours in a normal root canal.

pH



Results



Before treatment



After treatment



After 15 months



After 3 months



After 3 years



After 4 years



Aust Endod J 2021

ORIGINAL RESEARCH

Physicochemical properties and cytocompatibility of newly developed calcium silicate-based sealers

Min-Gyu Park, DDS¹; In-Ryoung Kim, PhD^{2,3}; Hyung Joon Kim, PhD^{3,4}; Sang Won Kwak, DDS, MS, PhD^{1,3}; and Hyeon-Cheol Kim, DDS, MS, PhD^{1,3}

¹ Department of Conservative Dentistry, Dental Research Institute, Pusan National University School of Dentistry, Yangsan, Korea

² Department of Oral Anatomy, Pusan National University School of Dentistry, Yangsan, Korea

³ Dental and Life Science Institute, Pusan National University, Yangsan, Korea

⁴ Department of Oral Physiology, Pusan National University School of Dentistry, Yangsan, Korea

⁵ Periodontal Diseases Signalling Network Research Center, Pusan National University, Yangsan, Korea



Cell Attachment Evaluation (One-Fil)



BIOCERAMIC ROOT CANAL SEALER FDA

CONTENTS

508-1 One-Fil 2g x 1ea

508-2 One-Fil 0.5g x 1ea