

Any-Com® Bulk

FDA
POSTERIOR BULK FILL FLOWABLE RESIN


It can be used as a base-liner for both direct and indirect restoration up to 4mm a single-layer placement. It is also suitable for lining deep cavity. (Class cavity)

Applications


- Base in class I, II direct restorations
- Liner under direct restorative materials
- Fluoride releasing resin
- Small class I
- Core build-up
- Post endodontic coronal seal
- Primary molars

Advantages


- Save time
- Safe filling
- Simple technology
- One layer curing up to 4mm
- Minimized polymerization shrinkage stress



Results



Core build-up Application






Resin Tag




Transparency Shade

Physical properties



Self-leveling



Time →



POSTERIOR BULK-FILL
FLOW RESIN FDA

CONTENTS

307-11~307-15 Any-Com Bulk 2g x 1ea



ETCHING

BONDING

RESTORATION