

Any-Com® Flow

FDA
NANO-HYBRID FLOWABLE COMPOSITE

Light-cured, radiopaque nano-hybrid flowable composite for the restorative treatment.

Applications

- Class III, IV, V restoration
- Base-liner
- Repair of composite resin, crown, bridge, and temporary restoration
- Pit & Fissure sealant

Advantages

- Excellent flowability and low shrinkage in cervical and proximal area
- Replaceable with sealant and base liner
- Easy to apply to narrow and deep carious areas
- Variety of shade choices
- Under-cut block out
- Depth of cure : 3mm (10 sec.)

Results



Clinical reference [1]



Clinical reference [2]




SHADE


A1	A2	A3	A3.5	A4
B1	B2	C2	OA2	OA3
OWT	TL			

Physical properties


Flexural strength (MPa)



Shrinkage rate (%)



Compressive strength (MPa)



NANO-HYBRID FLOWABLE COMPOSITE

CONTENTS

# 306H	Any-Com Flow Kit
	Any-Com Flow 2g x 4ea, Any-Etch 3ml x 1ea, Hi-Bond Universal 5ml x 1ea
# 306-00	Any-Com Flow Refill 2g x 2ea

ETCHING

BONDING

RESTORATION