

STRONG RETENTION EASY REMOVAL


HIGH RETENTION
EASY CLEAN-UP
EASY HANDLING

PRODUCT	WORKING TIME (at 25 °C)	SETTING TIME (at 25 °C)
Any-Temp NE	1'00''	3'00''
Any-Cem Implant	2'00''	3'10''


Advantages

- Dual-cure
- Excellent radiopacity
- Easy to remove the excess cement
- Excellent flowability
- Retrievable
- Low film thickness


Shear bond strength (MPa)



Compressive strength (MPa)



Radiopacity



Any-Cem Implant

For
long term
temporary



CONTENTS

- # 202-1 Any-Cem Implant Pink 6g x 1syringe
- # 202-2 Any-Cem Implant White 6g x 1syringe

Any-Temp NE

For
temporary



CONTENTS

- # 201-1 Any-Temp NE 10g x 2syringes
- # 201-2 Any-Temp NE 10g x 3syringes